



Strong. Aesthetic. Metal-free.

ZERAMEX®



+ swiss made

Ceramics – a confident choice!

ZERAMEX® is a pioneer among modern two-piece ceramic implants. As a passionate innovator, ZERAMEX® is constantly working on developments in implantology. In doing so, it builds on the Swiss tradition of processing hard zirconia, called the “white diamond,” from which the implant is produced. As a long-term partner, ZERAMEX® allows you a range of options and can enhance your practice through reliable products that are easy to use and are tailored to the needs. What this means:

- **Key expertise in ceramic processing and Swiss quality**
- **An established implant system with an impressive performance record**



10 years of experience
in two-piece ceramic
implants



ZERAMEX® life-long
guarantee of implants*



ZERAMEX® 10-year
guarantee of secondary
components*



Over 25,000 implants
sold



FDA approval



Focus on ceramic



Customer consultation by
ceramic experts



Innovative

* The current warranty conditions apply. These can be found on our website www.zeramex.com.

For maximum patient satisfaction.

The requirements placed upon dental implants have increased in recent years. First and foremost is the patient's wish to have both a healthy and attractive solution. ZERAMEX® is synonymous with a high quality of life thanks to the metal-free ceramic implants. The clinical use of ZERAMEX® implants can boast of an impressive success rate of over 96 percent⁹ healing. Ceramic's high resistance to corrosion combined with the low plaque affinity minimizes the risk of inflammation¹¹. In addition, blood flow to the gums around the implant is retained: One study found that the ceramic implant has characteristics similar to those of a natural tooth². The benefits of ZERAMEX®:

> Long-term aesthetics

> High compatibility



Metal-free and corrosion-resistant⁴:



Simply attractive¹



Minimize risks^{4,11,12,14}



Swiss made



Gums love ceramic^{2,11}



Clinically successful^{9,17}


The ZERAMEX® implant family.

ZERAMEX® OFFERS YOU THE IDEAL SOLUTION!

The new ZERAMEX® XT implant is the latest member in the family of two-piece, reversible screw-in ceramic implants. The successful ZERAMEX® P6 implant with a soft tissue level design is ideal in the posterior area and offers easy access to the implant. The root-shaped design of the ZERAMEX® XT implant achieves high primary stability, and high prosthetic flexibility is ensured thanks to the new internal connection.




ZERAMEX® XT



- Conical implant shape (root-shaped)
- Bolt-in tube internal connection
- Variable placement depth (optional up to 0.6 mm supracrestal)
- ATZ ceramics
- VICARBO® screw




ZERAMEX® P6




- Tissue Level Design
- Hexagon head
- Variable placement depth (optional up to 0.8 mm supracrestal)
- ATZ ceramics
- VICARBO® screw

Strong.


The latest generation of ZERAMEX® ceramic implants achieves strong connections thanks to the carbon-ceramic technology. The key component is the VICARBO® screw made of carbon-fiber-reinforced high-tech plastic. This material has excellent isoelastic properties and fits tightly against the existing contour. It creates a high-strength and form-fitting, albeit tension-free screwed ceramic-ceramic connection. The design of the connection was developed with the aim of minimizing bacterial colonization and accumulation.



The VICARBO® screw precisely fits the existing contour of the thread upon tightening, thanks to the very different E-modulus.



The continuous, form-fitting carbon fibers consistently absorb tensile forces and give the VICARBO® screw its high tensile strength.



VICARBO® HIGH-PERFORMANCE CARBON TECHNOLOGY


The carbon-fiber-reinforced plastic has excellent isoelastic properties and obvious strengths: Metal-free, biocompatible in compliance with ISO10993 for long-term implants. Modulus of elasticity: >160 GPa., tensile strength: 2000 MPa.¹⁵ (compared with Grade 5 titanium: Modulus of elasticity: 114 GPa., tensile strength: 1100 MPa.¹⁶)

OUTSTANDING MATERIAL STRENGTH

ZERAMEX® ceramic implants have a high fatigue strength. Implants and abutments are made from highly compressed, hot isostatic post-compacted (HIP) zirconium dioxide ATZ blanks (hard zirconia). The result is high strength, biomechanical stability and a precise fit.

No thermal process (sintering) and no finishing occur after the final shaping of the external and internal implant geometry. This ensures a high level of precision, and no more changes can occur in the material structure. This manufacturing process is very elaborate and requires much experience and know-how.

FATIGUE STRENGTHS IN COMPARISON (according to ISO 14801)¹⁷




A ZERAMEX® P6, ø 3.3 mm, SN, and ZERAMEX® P6 Abutment SN Straight

B ZERAMEX® P6, ø 3.3 mm, SN, and ZERAMEX® P6 Abutment SN Angular

Aesthetic.

The patient demand for ceramic implant solutions is continuously increasing. The megatrend towards white, metal-free dentistry, which started years ago, continues to grow. The two-piece, reversible screw connection can facilitate unhindered healing of the implant and optimize accommodation by soft tissue. Together with the white color and favorable tolerability by the gums¹¹, long-lasting aesthetics are ensured.




OPTIMUM GINGIVAL BLOOD FLOW

Blood circulation in the gingiva around zirconium dioxide is comparable with that around a natural tooth, while blood circulation in the gingiva around titanium is significantly lower.² The reduced build-up of plaque makes hygiene easier for patients and promotes inflammation-free gums.¹¹

NATURAL AESTHETICS

The white tooth color of zirconium dioxide is visually superior to gray titanium because there are no gray edges, and a dark implant core does not shine through.¹ In the case of thin gingiva or gingival retraction, ceramic implants have an advantage with their white tooth aesthetics.



Situation after potential recession of gingiva.

Left: Restoration with ZERAMEX®.

Right: restoration with a conventional titanium implant.

PROSTHETIC FLEXIBILITY

Thanks to the two-piece, reversible screw-in design, ZERAMEX® implants have high prosthetic flexibility, yet remain simple and well-designed systems. Use common procedures and provide full-ceramic, screw-in restorations. Through familiarization with the CAD/CAM systems of Exocad and 3Shape (ZERAMEX® P6), you will also benefit from the advantages of a digital workflow.

Metal-free.


The 100-percent metal-free ceramic implants closely approximate the natural tooth root in terms of aesthetics and function. They are well tolerated, are completely free of metallic corrosion, and conduct neither electricity nor heat. Titanium can release titanium ions, which accumulate in the surrounding soft tissue and can cause inflammation.⁵ Likewise, titanium's resistance to corrosion can decrease due to surface bacteria.¹⁴ The use of metal-free ceramic implants precludes these effects.

OUTSTANDING OSSEOUS INTEGRATION

The surface structure is optimized through sandblasting and etching to permit osteoblasts to grow directly on the implant, allowing a solid attachment to the surface of the implant.¹³



Sectional view after eight weeks: successful osseous integration of the hydrophilic ZERAFIL™ surface to ceramic. Working group of Prof. D. Buser, University of Bern, Switzerland¹³



MINIMIZING RISK FACTORS

Long-term studies demonstrate that peri-implantitis is a risk associated with implant treatments.³ If left untreated, peri-implantitis can result in the loss of the implant. The ZERAMEX® P6 implant can help minimize certain risk factors of peri-implantitis:

Plaque¹²:

→ Low plaque affinity of ceramic¹¹

Metallic corrosion^{4,14}:

→ Ceramic cannot succumb to metallic corrosion

Our recommendation:

IMPLANT TYPE	ZERAMEX® P6	ZERAMEX® XT
Anterior region	+	++
Posture region	++	+
Thin mucosal biotype	++	++
Vertical/lateral bone growth with simultaneous implant insertion (at least 5 mm residual bone height)	+	++
External sinus lift	possible	++
Covered healing	+	++
Open healing	++	+
Immediate implant insertion in extraction alveoli	+	++

+ recommended ++ highly recommended

You can use ZERAMEX® to treat the following indications:



Single tooth



Partial bridge



Full bridge

Order hotline DE/CH/AT

T 00800 93 55 66 37

F 00800 93 55 63 77

order@zeramex.com

Head office

Dentalpoint AG

Bodenackerstrasse 5

CH-8048 Zurich, Switzerland

T 0041 44 388 36 36

F 0041 44 388 36 39

European branch

Dentalpoint Germany GmbH

Wallbrunnstrasse 24

D-79539 Lörrach, Germany

T 0049 7621 1612 749

F 0049 7621 1612 780

info@zeramex.com

www.zeramex.com

References

- 1) Cosgarea R et al., Peri-implant soft tissue colour around titanium and zirconia abutments: a prospective randomized controlled clinical study. *Clinical Oral Implant Research* 26, 2015 / 537–544.
- 2) Kajiwaru N et al., Soft tissue biological response to zirconia and metal implant abutments compared with natural tooth: Microcirculation Monitoring as a Novel Bioindicator., *Implant Dentistry Volume* 24, Number 1 2015.
- 3) Derks J et al., Effectiveness of Implant Therapy Analyzed in a Swedish Population: Prevalence of Peri-implantitis. *J Dent Res.* 2016 Jan; 95(1):43-9.
- 4) Wachi T et al., Release of titanium ions from an implant surface and their effect on cytokine production related to alveolar bone resorption. *Toxicology.* 2015 Jan 2; 327:1-9.
- 5) Addison O et al., Do 'passive' medical titanium surfaces deteriorate in service in the absence of wear? *J R Soc Interface.* 2012 Nov 7; 9(76):3161-4.
- 6) Derks J et al., Effectiveness of implant therapy analyzed in a Swedish population: early and late implant loss. *J Dent Res.* 2015 Mar; 94(3 Suppl):44S-51S.
- 7) Oliva J et al., Five-year success rate of 831 consecutively placed Zirconia dental implants in humans: a comparison of three different rough surfaces. *Int J Oral Maxillofac Implants.* 2010 Mar-Apr;25(2):336-44.
- 8) Manzano G et al., Comparison of clinical performance of zirconia implants and titanium implants in animal models: a systematic review. *Int J Oral Maxillofac Implants.* 2014 Mar-Apr;29(2):311-20.
- 9) Jank S et al., Success Rate of Two-Piece Zirconia Implants: A Retrospective Statistical Analysis. *Implant Dent.* 2016 Feb 1.
- 10) Cionca N, Two-piece zirconia implants supporting all-ceramic crowns: a prospective clinical study. *Clin Oral Implants Res.* 2015 Apr;26(4):413-8.
- 11) Scarano A et al., Bacterial adhesion on commercially pure titanium and zirconium oxide disks: an in vivo human study. *J Periodontol.* 2004 Feb; 75(2):292-6.
- 12) Canullo L et al., Distinguishing predictive profiles for patient-based risk assessment and diagnostics of plaque-induced, surgically and prosthetically triggered peri-implantitis. *Clin Oral Implants Res.* 2015 Nov 20.
- 13) Chappuis V et al., Osseointegration of zirconia and titanium implants in the presence of multinucleated giant cells. *CIDRR*, 2015 Sept. 17.
- 14) Sridhar S et al., In Vitro Investigation of the Effect of Oral Bacteria in the Surface Oxidation of Dental Implants. *Clin Implant Dent Relat Res.* 2015 Oct;17 Suppl 2:e562-75.
- 15) Invisio T-PB-PUR-E-0031-A (3/2013).
- 16) Boyer R et al., *Materials Properties Handbook: Titanium Alloys*, ASM International, 1994.
- 17) Ermüdungstests nach ISO14801 [Fatigue test in compliance with IOS14801]; Report No. 16010106-D-CS of 03/31/2016 and Report No. 14070102-D-CS of 04/21/2015; Study director: Nicolas Graf; Spineserv GmbH & Co. KG, Söflinger Strasse 100, D-89077 Ulm

© 2017 Dentalpoint AG—All rights reserved. ZERAMEX®, ZERAFIX™, ZERALOCK™ and VICARBO® are trademarks or registered trademarks of Dentalpoint AG.

Swiss Dental Solutions® is a registered trademark of Swiss Dental Solutions AG.

Locator® is a registered trademark of Zest Anchors, Inc.

Exocad is a registered trademark of Exocad GmbH, Germany.

3Shape is a registered trademark of 3Shape A/S, Denmark.

ZERAMEX®

Strong. Aesthetic. Metal-free.