

In rapid succession: Zeramex launches further innovations

No compromises, natural aesthetics and the highest prosthetic flexibility: under this motto, Zeramex is launching the Zeramex XT 3.5 mm ceramic implant, the first reversibly screwable and 100 percent metal-free abutments for removable dentures and the new, smaller Zerabase X adhesive base.

1a Zeramex XT
Implant systems
Implants with an
endosseous diameter
of 3.5 mm



1 Small-Base Implants



2 Individual Abutment



3 to 4 Standard-Abutments


The portfolio of the 100 percent metal-free and reversibly screwable Zeramex XT ceramic implant system is now complete. After the launch of the new competence center Zeramex Digital Solutions, which serves the demand for individual abutments and one-piece, monolithic crowns with implant connection, Firrna Dentalpoint is launching further innovations on the market.

The portfolio of the Zeramex XT implant system now also includes the anterior teeth in the lower jaw and lateral incisors in the upper jaw

Implants with an endosseous diameter of 3.5 mm (Fig. 1a). The new small-base implants (SB) (Fig. 1) are available in lengths of 8, 10 and 12 mm. In addition to the familiar regular and wide-base platforms with endosseous diameters of 4.2 mm for the RB implant and 5.5 mm for the WB implant, the range of Zeramex XT implants is now complete (Figs. 2 to 5).

In addition to the SB implants, the system has also been expanded with new abutments for removable dentures. The Zeramex Docklocs abutments are the first 100

5 Product family Zeramex



6 to 9 Zeramex Docklocs Abutments

Percent metal-free and reversibly screwable abutments for removable dentures. They are available in heights of 2, 3 and 4 mm and fit on all Zeramex XT SB, RB and WB implants (Fig 6 to 9).



In order to further optimize the digital workflow, a smaller adhesive base has recently become available. Like the Zerabase adhesive base, the new Zerabase X is also available with and without abutment feet (Figs. 10 and 11) for single crowns or bridge restorations (engaged / unengaged) as well as for all SB, RB and WB platforms.

"Zeramex now offers the most comprehensive range in ceramic implantology to achieve natural aesthetic results," explains Dr. Ricarda Jansen, Director of Dental Implants at Zeramex.

The innovations SB implants, Zeramex Docklocs abutments and the new Zerabase X adhesive base are available now.

Dentalpoint AG / Zeramex
 BodenackerstraBe 5
 8957 Spreitenbach / Switzerland
 Phone 00800 93556637 Mail
 info@zeramex.com

10 Zerabase X with Abutment-foot



11 Zerabase X one Abutment tube